



TNFRSF10A Antibody, Biotin conjugated

Product Code	CSB-PA023964LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00220
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 10A protein (263-468AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 10A (Death receptor 4) (TNF-related apoptosis-inducing ligand receptor 1) (TRAIL receptor 1) (TRAIL-R1) (CD antigen CD261), TNFRSF10A, APO2 DR4 TRAILR1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TNFRSF10A