



TNMD Antibody, FITC conjugated

Product Code	CSB-PA024007LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2S6
Immunogen	Recombinant Human Tenomodulin protein (131-254AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tenomodulin (TeM) (hTeM) (Chondromodulin-1-like protein) (ChM1L) (hChM1L) (Chondromodulin-I-like protein) (Myodulin) (Tendin), TNMD, CHM1L
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TNMD