



TRIM21 Antibody, Biotin conjugated

Product Code	CSB-PA024457LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19474
Immunogen	Recombinant Human E3 ubiquitin-protein ligase TRIM21 protein (1-475AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase TRIM21 (EC 2.3.2.27) (52 kDa Ro protein) (52 kDa ribonucleoprotein autoantigen Ro/SS-A) (RING finger protein 81) (RING-type E3 ubiquitin transferase TRIM21) (Ro(SS-A)) (Sjogren syndrome type A antigen) (SS-A) (Tripartite motif-containing protein 21), TRIM21, RNF81 RO52 SSA1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TRIM21
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.