





TRPC1 Antibody

Product Code	CSB-PA025075LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48995
Immunogen	Recombinant Human Short transient receptor potential channel 1 protein (39-207AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Short transient receptor potential channel 1 (TrpC1) (Transient receptor protein 1) (TRP-1), TRPC1, TRP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TRPC1

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA025075LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA025075LA01HU at dilution of 1:100

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.