



IKBKB Antibody, Biotin conjugated

Product Code	CSB-PA16309D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14920
Immunogen	Recombinant Human Inhibitor of nuclear factor kappa-B kinase subunit beta protein (395-588AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inhibitor of nuclear factor kappa-B kinase subunit beta (I-kappa-B-kinase beta) (IKK-B) (IKK-beta) (IkBKB) (EC 2.7.11.10) (I-kappa-B kinase 2) (IKK2) (Nuclear factor NF-kappa-B inhibitor kinase beta) (NFKBIBK), IKBKB, IKKB
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	IKBKB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.