



tal1 Antibody, Biotin conjugated

Product Code	CSB-PA20969D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O93507
Immunogen	Recombinant Zebrafish T-cell acute lymphocytic leukemia protein 1 homolog protein (1-324AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell acute lymphocytic leukemia protein 1 homolog (TAL-1) (Stem cell protein) (zSCL), tal1, scl tal-1
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Others
Target Names	tal1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.