



NP Antibody, Biotin conjugated

Product Code	CSB-PA325769LD01IIM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15682
Immunogen	Recombinant Influenza A virus Nucleoprotein protein (1-498AA)
Raised In	Rabbit
Species Reactivity	Influenza A virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoprotein (Nucleocapsid protein) (Protein N), NP
Immunogen Species	Influenza A virus (strain A/Wilson-Smith/1933 H1N1) (Influenza A virus (strain A/WS/1933 H1N1))
Research Area	Others
Target Names	NP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.