




# Thiaminase-1 Antibody

<b>Product Code</b>	CSB-PA332698HA01PAI
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P45741
<b>Immunogen</b>	Recombinant Paenibacillus thiaminolyticus Thiaminase-1 protein (30-409AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Paenibacillus thiaminolyticus
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Thiaminase-1 (EC 2.5.1.2) (Thiaminase I) (Thiamine pyridinylase)
<b>Immunogen Species</b>	Paenibacillus thiaminolyticus
<b>Research Area</b>	Others

## Image



### Western Blot

Positive WB detected in Recombinant protein

All lanes: Thiaminase-1 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 58 kDa

Observed band size: 58 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.