



Thiaminase-1 Antibody, Biotin conjugated

Product Code	CSB-PA332698HD01PAI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45741
Immunogen	Recombinant Paenibacillus thiaminolyticus Thiaminase-1 protein (30-409AA)
Raised In	Rabbit
Species Reactivity	Paenibacillus thiaminolyticus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thiaminase-1 (EC 2.5.1.2) (Thiaminase I) (Thiamine pyridinylase)
Immunogen Species	Paenibacillus thiaminolyticus (Bacillus thiaminolyticus)
Research Area	Others