




ADK Antibody

Product Code	CSB-PA346232LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55263
Immunogen	Recombinant Human Adenosine kinase protein (18-256AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adenosine kinase (AK) (EC 2.7.1.20) (Adenosine 5'-phosphotransferase), ADK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADK

Image



Western Blot

Positive WB detected in: HL60 whole cell lysate, Rat kidney tissue, Rat liver tissue, Mouse kidney tissue, Mouse liver tissue


All lanes: ADK antibody at 3µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 41, 39, 35, 37 kDa


Observed band size: 41 kDa



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA346232LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA346232LA01HU at dilution of 1:100



Immunofluorescent analysis of PC-3 cells using CSB-PA346232LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.