




btfP Antibody

Product Code	CSB-PA346537LA01BDP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54355
Immunogen	Recombinant Bacteroides fragilis Fragilysin protein (26-405AA)
Raised In	Rabbit
Species Reactivity	Bacteroides fragilis
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fragilysin (EC 3.4.24.74) (Enterotoxin), btfP
Immunogen Species	Bacteroides fragilis
Research Area	Others
Target Names	btfP

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: btfP antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 58 kDa

Observed band size: 58 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.