



btFP Antibody, Biotin conjugated

Product Code	CSB-PA346537LD01BDP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54355
Immunogen	Recombinant Bacteroides fragilis Fragilysin protein (26-405AA)
Raised In	Rabbit
Species Reactivity	Bacteroides fragilis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fragilysin (EC 3.4.24.74) (Enterotoxin), btFP
Immunogen Species	Bacteroides fragilis
Research Area	Others
Target Names	btFP