



btFP Antibody, Biotin conjugated

Product Code	CSB-PA346537LD01BDP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54355
Immunogen	Recombinant Bacteroides fragilis Fragilysin protein (26-405AA)
Raised In	Rabbit
Species Reactivity	Bacteroides fragilis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fragilysin (EC 3.4.24.74) (Enterotoxin), btFP
Immunogen Species	Bacteroides fragilis
Research Area	Others
Target Names	btFP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.