



phr Antibody, HRP conjugated

Product Code	CSB-PA361494LB01FPZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05327
Immunogen	Recombinant Synechococcus sp. Deoxyribodipyrimidine photo-lyase protein (2-484AA)
Raised In	Rabbit
Species Reactivity	Synechococcus sp.
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deoxyribodipyrimidine photo-lyase (EC 4.1.99.3) (DNA photolyase) (Photoreactivating enzyme), phr, phrA
Immunogen Species	Synechococcus sp. (strain ATCC 27144 / PCC 6301 / SAUG 1402/1) (Anacystis nidulans)
Research Area	Others
Target Names	phr
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.