




THI2.1 Antibody

Product Code	CSB-PA365586LA01VDO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01538
Immunogen	Recombinant Viscum album Viscotoxin-A3 protein (27-72AA)
Raised In	Rabbit
Species Reactivity	Viscum album
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Viscotoxin-A3 [Cleaved into: Viscotoxin-A3; Acidic protein], THI2.1
Immunogen Species	Viscum album
Research Area	Others
Target Names	THI2.1

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: THI2.1 antibody at 2.7μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 20 kDa

Observed band size: 20 kDa