



kgp Antibody, Biotin conjugated

Product Code	CSB-PA464342LD01EXZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B2RLK2
Immunogen	Recombinant Porphyromonas gingivalis Lys-gingipain protein (229-594AA)
Raised In	Rabbit
Species Reactivity	Porphyromonas gingivalis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lys-gingipain (EC 3.4.22.47) (Lysine-specific cysteine proteinase Kgp) [Cleaved into: Lys-gingipain catalytic subunit; 39 kDa adhesin; 15 kDa adhesin; 44 kDa adhesin], kgp
Immunogen Species	Porphyromonas gingivalis (strain ATCC 33277 / DSM 20709 / CIP 103683 / JCM 12257 / NCTC 11834 / 2561)
Research Area	Others
Target Names	kgp
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.