



ghrA Antibody, Biotin conjugated

Product Code	CSB-PA482379HD01ENM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B7LFE3
Immunogen	Recombinant Escherichia coli Glyoxylate/hydroxypyruvate reductase A protein (1-312AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glyoxylate/hydroxypyruvate reductase A (EC 1.1.1.79) (EC 1.1.1.81) (2-ketoacid reductase), ghrA
Immunogen Species	Escherichia coli (strain 55989 / EAEC)
Research Area	Others
Target Names	ghrA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.