




Vertebrate ancient opsin Antibody

Product Code	CSB-PA517517LA01SWI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O13018
Immunogen	Recombinant Salmo salar Vertebrate ancient opsin protein (1-75AA)
Raised In	Rabbit
Species Reactivity	Salmo salar
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vertebrate ancient opsin
Immunogen Species	Salmo salar
Research Area	Others

Image



Western Blot

Positive WB detected in Recombinant protein
All lanes: Vertebrate ancient opsin antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 25 kDa

Observed band size: 25 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.