



# Vertebrate ancient opsin Antibody, FITC conjugated

Product Code	CSB-PA517517LC01SWI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O13018
Immunogen	Recombinant Salmo salar Vertebrate ancient opsin protein (1-75AA)
Raised In	Rabbit
Species Reactivity	Salmo salar
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vertebrate ancient opsin
Immunogen Species	Salmo salar (Atlantic salmon)
Research Area	Others