



Vertebrate ancient opsin Antibody, Biotin conjugated

Product Code	CSB-PA517517LD01SWI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O13018
Immunogen	Recombinant Salmo salar Vertebrate ancient opsin protein (1-75AA)
Raised In	Rabbit
Species Reactivity	Salmo salar
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vertebrate ancient opsin
Immunogen Species	Salmo salar (Atlantic salmon)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.