



# KDM6B Antibody, Biotin conjugated

Product Code	CSB-PA517674LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15054
Immunogen	Recombinant Human Lysine-specific demethylase 6B protein (1530-1682AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysine-specific demethylase 6B (EC 1.14.11.-) (JmjC domain-containing protein 3) (Jumonji domain-containing protein 3) (Lysine demethylase 6B), KDM6B, JMJD3 KIAA0346
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	KDM6B