



KDM6B Antibody, Biotin conjugated

Product Code	CSB-PA517674LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15054
Immunogen	Recombinant Human Lysine-specific demethylase 6B protein (1530-1682AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysine-specific demethylase 6B (EC 1.14.11.-) (JmjC domain-containing protein 3) (Jumonji domain-containing protein 3) (Lysine demethylase 6B), KDM6B, JMJD3 KIAA0346
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	KDM6B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.