



PTK2B Antibody, Biotin conjugated

Product Code	CSB-PA613494LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14289
Immunogen	Recombinant Human Protein-tyrosine kinase 2-beta protein (707-875AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein-tyrosine kinase 2-beta (EC 2.7.10.2) (Calcium-dependent tyrosine kinase) (CADTK) (Calcium-regulated non-receptor proline-rich tyrosine kinase) (Cell adhesion kinase beta) (CAK-beta) (CAKB) (Focal adhesion kinase 2) (FADK 2) (Proline-rich tyrosine kinase 2) (Related adhesion focal tyrosine kinase) (RAFTK), PTK2B, FAK2 PYK2 RAFTK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTK2B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.