



PDK1 Antibody, Biotin conjugated

Product Code	CSB-PA613584LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15118
Immunogen	Recombinant Human [Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 1, mitochondrial protein (294-432AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	[Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 1, mitochondrial (EC 2.7.11.2) (Pyruvate dehydrogenase kinase isoform 1) (PDH kinase 1), PDK1, PDHK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.