



LDB1 Antibody, HRP conjugated

Product Code	CSB-PA771433LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86U70
Immunogen	Recombinant Human LIM domain-binding protein 1 protein (50-375AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LIM domain-binding protein 1 (LDB-1) (Carboxyl-terminal LIM domain-binding protein 2) (CLIM-2) (LIM domain-binding factor CLIM2) (hLdb1) (Nuclear LIM interactor), LDB1, CLIM2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LDB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.