



SGMS1 Antibody, Biotin conjugated

Product Code	CSB-PA801243LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86VZ5
Immunogen	Recombinant Human Phosphatidylcholine:ceramide cholinephosphotransferase 1 protein (48-137AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylcholine:ceramide cholinephosphotransferase 1 (EC 2.7.8.27) (Medulla oblongata-derived protein) (Protein Mob) (Sphingomyelin synthase 1) (Transmembrane protein 23), SGMS1, MOB SMS1 TMEM23
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SGMS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.