



UHK1 Antibody, Biotin conjugated

Product Code	CSB-PA840554LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAS1
Immunogen	Recombinant Human Serine/threonine-protein kinase Kist protein (183-385AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase Kist (EC 2.7.11.1) (Kinase interacting with stathmin) (PAM COOH-terminal interactor protein 2) (P-CIP2) (U2AF homology motif kinase 1), UHK1, KIS KIST
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UHK1