



Fuca1 Antibody, HRP conjugated

| | |
|---------------------|--|
| Product Code | CSB-PA858759LB01MO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q99LJ1 |
| Immunogen | Recombinant Mouse Tissue alpha-L-fucosidase protein (18-452AA) |
| Raised In | Rabbit |
| Species Reactivity | Mouse |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Tissue alpha-L-fucosidase (EC 3.2.1.51) (Alpha-L-fucosidase I) (Alpha-L-fucoside fucosylhydrolase 1) (Alpha-L-fucosidase 1), Fuca1, Fuca |
| Immunogen Species | Mus musculus (Mouse) |
| Research Area | Signal Transduction |
| Target Names | Fuca1 |