



Fuca1 Antibody, HRP conjugated

Product Code	CSB-PA858759LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99LJ1
Immunogen	Recombinant Mouse Tissue alpha-L-fucosidase protein (18-452AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tissue alpha-L-fucosidase (EC 3.2.1.51) (Alpha-L-fucosidase I) (Alpha-L-fucoside fucohydrolase 1) (Alpha-L-fucosidase 1), Fuca1, Fuca
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Fuca1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.