



EPAS1 Antibody, HRP conjugated

Product Code	CSB-PA860780LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99814
Immunogen	Recombinant Human Endothelial PAS domain-containing protein 1 protein (543-870AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endothelial PAS domain-containing protein 1 (EPAS-1) (Basic-helix-loop-helix-PAS protein MOP2) (Class E basic helix-loop-helix protein 73) (bHLHe73) (HIF-1-alpha-like factor) (HLF) (Hypoxia-inducible factor 2-alpha) (HIF-2-alpha) (HIF2-alpha) (Member of PAS protein 2) (PAS domain-containing protein 2), EPAS1, BHLHE73 HIF2A MOP2 PASD2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	EPAS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.