



SLC32A1 Antibody, Biotin conjugated

Product Code	CSB-PA863984LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H598
Immunogen	Recombinant Human Vesicular inhibitory amino acid transporter protein (14-118AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vesicular inhibitory amino acid transporter (GABA and glycine transporter) (Solute carrier family 32 member 1) (Vesicular GABA transporter) (hVIAAT), SLC32A1, VGAT VIAAT
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC32A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.