





DIABLO Antibody

Product Code	CSB-PA865113LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NR28
Immunogen	Recombinant Human Diablo homolog, mitochondrial protein (168-236AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diablo homolog, mitochondrial (Direct IAP-binding protein with low pI) (Second mitochondria-derived activator of caspase) (Smac), DIABLO, SMAC
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DIABLO

Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA865113LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA865113LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)