




Pectate lyase 1 Antibody

Product Code	CSB-PA871029LA01COAD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9SCG9
Immunogen	Recombinant Hesperocyparis arizonica Pectate lyase 1 protein (22-367AA)
Raised In	Rabbit
Species Reactivity	Hesperocyparis arizonica
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pectate lyase 1 (EC 4.2.2.2) (Major pollen allergen Cup a 1) (allergen Cup a 1)
Immunogen Species	Hesperocyparis arizonica
Research Area	Others

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: Pectate lyase 1 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 43 kDa

Observed band size: 43 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.