



Pectate lyase 1 Antibody, Biotin conjugated

Product Code	CSB-PA871029LD01COAD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9SCG9
Immunogen	Recombinant Hesperocyparis arizonica Pectate lyase 1 protein (22-367AA)
Raised In	Rabbit
Species Reactivity	Hesperocyparis arizonica
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pectate lyase 1 (EC 4.2.2.2) (Major pollen allergen Cup a 1) (allergen Cup a 1)
Immunogen Species	Hesperocyparis arizonica (Arizona cypress) (Cupressus arizonica)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.