



HIPK3 Antibody, HRP conjugated

Product Code	CSB-PA872492LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H422
Immunogen	Recombinant Human Homeodomain-interacting protein kinase 3 protein (841-1130AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeodomain-interacting protein kinase 3 (EC 2.7.11.1) (Androgen receptor-interacting nuclear protein kinase) (ANPK) (Fas-interacting serine/threonine-protein kinase) (FIST) (Homolog of protein kinase YAK1), HIPK3, DYRK6 FIST3 PKY
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	HIPK3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.