



SENP2 Antibody, Biotin conjugated

Product Code	CSB-PA888026LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HC62
Immunogen	Recombinant Human Sentrin-specific protease 2 protein (138-324AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sentrin-specific protease 2 (EC 3.4.22.-) (Axam2) (SMT3-specific isopeptidase 2) (Smt3ip2) (Sentrin/SUMO-specific protease SENP2), SENP2, KIAA1331
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SENP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.