



Snaclec rhodocytin subunit alpha Antibody, HRP conjugated

Product Code	CSB-PA888242LB01CBG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9I841
Immunogen	Recombinant Calloselasma rhodostoma Snaclec rhodocytin subunit alpha protein (1-136AA)
Raised In	Rabbit
Species Reactivity	Calloselasma rhodostoma
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Snaclec rhodocytin subunit alpha (Aggretin alpha chain) (Rhodoaggretin subunit alpha)
Immunogen Species	Calloselasma rhodostoma (Malayan pit viper) (Agkistrodon rhodostoma)
Research Area	Others