



Snaclec rhodocytin subunit alpha Antibody, FITC conjugated

Product Code	CSB-PA888242LC01CBG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9I841
Immunogen	Recombinant Calloselasma rhodostoma Snaclec rhodocytin subunit alpha protein (1-136AA)
Raised In	Rabbit
Species Reactivity	Calloselasma rhodostoma
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Snaclec rhodocytin subunit alpha (Aggretin alpha chain) (Rhodoaggretin subunit alpha)
Immunogen Species	Calloselasma rhodostoma (Malayan pit viper) (Agkistrodon rhodostoma)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.