



ABCG2 Antibody, FITC conjugated

Product Code	CSB-PA891568LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNQ0
Immunogen	Recombinant Human ATP-binding cassette sub-family G member 2 protein (296-371AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-binding cassette sub-family G member 2 (Breast cancer resistance protein) (CDw338) (Mitoxantrone resistance-associated protein) (Placenta-specific ATP-binding cassette transporter) (Urate exporter) (CD antigen CD338), ABCG2, ABCP BCRP BCRP1 MXR
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ABCG2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.