





BCL2L11 Antibody

Product Code	CSB-PA002615LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43521
Immunogen	Recombinant Human Bcl-2-like protein 11 protein (1-198AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bcl-2-like protein 11 (Bcl2-L-11) (Bcl2-interacting mediator of cell death), BCL2L11, BIM
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BCL2L11


Image



Immunohistochemistry analysis of human testis tissue using CSB-PA002615LA01HU at dilution of 1:100



Immunocytochemistry analysis of human lung cancer using CSB-PA002615LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA002615LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.