



BCL2L11 Antibody, Biotin conjugated

Product Code	CSB-PA002615LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43521
Immunogen	Recombinant Human Bcl-2-like protein 11 protein (1-198AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bcl-2-like protein 11 (Bcl2-L-11) (Bcl2-interacting mediator of cell death), BCL2L11, BIM
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BCL2L11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.