





YBEY Antibody

Product Code	CSB-PA003773LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P58557
Immunogen	Recombinant Human Endoribonuclease YbeY protein (54-167AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoribonuclease YbeY (EC 3.1.-.-), YBEY, C21orf57
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	YBEY

Image



Immunohistochemistry analysis of human colon cancer using CSB-PA003773LA01HU at dilution of 1:100



Immunocytochemistry analysis of human tonsil tissue using CSB-PA003773LA01HU at dilution of 1:100