



CDH15 Antibody, Biotin conjugated

Product Code	CSB-PA005040LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55291
Immunogen	Recombinant Human Cadherin-15 protein (639-766AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cadherin-15 (Cadherin-14) (Muscle cadherin) (M-cadherin), CDH15, CDH14 CDH3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CDH15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.