



EGFR Antibody, HRP conjugated

Product Code	CSB-PA005571LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00533
Immunogen	Recombinant Human Epidermal growth factor receptor protein (1067-1210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Epidermal growth factor receptor (EC 2.7.10.1) (Proto-oncogene c-ErbB-1) (Receptor tyrosine-protein kinase erbB-1), EGFR, ERBB ERBB1 HER1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EGFR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.