




gfap Antibody

Product Code	CSB-PA009369LA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q58EE9
Immunogen	Recombinant Zebrafish Glial fibrillary acidic protein (1-444AA)
Raised In	Rabbit
Species Reactivity	Zebrafish, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glial fibrillary acidic protein (GFAP), gfap
Immunogen Species	Zebrafish
Research Area	Neuroscience
Target Names	gfap

Image



Western Blot

Positive WB detected in: MCF-7 whole cell lysate

All lanes: GFAP antibody at 3.4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 50, 51 kDa

Observed band size: 50 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.