



gfap Antibody, FITC conjugated

Product Code	CSB-PA009369LC01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q58EE9
Immunogen	Recombinant Zebrafish Glial fibrillary acidic protein (1-444AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glial fibrillary acidic protein (GFAP), gfap
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Neuroscience
Target Names	gfap
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.