



HTR1A Antibody, Biotin conjugated

Product Code	CSB-PA010881LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08908
Immunogen	Recombinant Human 5-hydroxytryptamine receptor 1A protein (230-334AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	5-hydroxytryptamine receptor 1A (5-HT-1A) (5-HT1A) (G-21) (Serotonin receptor 1A), HTR1A, ADRB2RL1 ADRBRL1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HTR1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.