



MYF5 Antibody, HRP conjugated

Product Code	CSB-PA015287LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13349
Immunogen	Recombinant Human Myogenic factor 5 protein (1-255AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myogenic factor 5 (Myf-5) (Class C basic helix-loop-helix protein 2) (bHLHc2), MYF5, BHLHC2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MYF5