



NCF2 Antibody, HRP conjugated

Product Code	CSB-PA015528LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19878
Immunogen	Recombinant Human Neutrophil cytosol factor 2 protein (178-526AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neutrophil cytosol factor 2 (NCF-2) (67 kDa neutrophil oxidase factor) (NADPH oxidase activator 2) (Neutrophil NADPH oxidase factor 2) (p67-phox), NCF2, NOXA2 P67PHOX
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NCF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.