



NCF2 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA015528LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P19878 |
| Immunogen | Recombinant Human Neutrophil cytosol factor 2 protein (178-526AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Neutrophil cytosol factor 2 (NCF-2) (67 kDa neutrophil oxidase factor) (NADPH oxidase activator 2) (Neutrophil NADPH oxidase factor 2) (p67-phox), NCF2, NOXA2 P67PHOX |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | NCF2 |